


Productivity Series Spectrophotometer

Model: S-1200E

- UNICO 120 0SERIES VISIBLE SPECTROPHOTOMETER
- 325-1000 NM WAVELENGTH RANGE
- 5NM SPECTRAL BANDWIDTH
- SAMPLE HOLDERS FOR 10 MM SINGLE SQUARE CUVETTE AND V-TYPE HOLDER FOR ROUND TUBE
- USB PORT FOR OPTIONAL SOFTWARE
- INCLUDED ARE 2PCS OF 10MM SQUARE GLASS CUVETTES; 12 PCS OF 10 MM TEST TUBES; DUST COVER; USER MANUAL
- VOLTAGE PRESET AT 220V



UNICO 1200 Series are reliable and expandable. UNICO Model 1200 has Absorbance, Transmittance and Concentration modes with automatic zeroing and blanking in a single button push. Its optional features include Window based application software for easy data collection by using Microsoft Excel. Its tool-free and alignment-free lamp is ideal for any standard spectrophotometer application, while staying within a conservative budget.



ADDITIONAL INFORMATION

<p>ACCESSORIES</p>	<p>Software</p> <p>S-1200-401 Unico S-1200 Spectrophotometer Software. Win 7 Or Win 10 Required. Functions include Standard Curve, Abs./%T/Conc. and Abs. vs. Time. includes USB Cable.</p> <p>S-1205-401 Unico S-1205 Spectrophotometer Software. Win 7 Or Win 10 Required. Functions include Standard Curve, Abs./%T/Conc., and Kinetics and Scanning. includes USB Cable.</p> <p>Accessories</p> <p>S-1200-101 Test Tube V-Type Cell Holder with Base for 8~25 Mm Diameter Test Tubes for S-1200 Series</p> <p>S-1200-102 Single Square Cuvette Holder with Base for 10 Mm Cuvettes for S-1200 Series</p> <p>S-1200-103 4-Position Cell Holder for 10 Mm Square Cuvettes (Standard in Model 1201) for S-1200 Series</p> <p>S-1200-104 Long Path 4-Cell Holder for Up To 50 Mm Rectangular Cuvettes for S-1200 Series</p> <p>S-1200-105 Water-Jacketed Cell Holder W/Base and Panel for Temperature Controlled Application for S-1200 Series</p> <p>S-1200-106 Universal Base for V-Type Holder, Single Cell Holder and Water-Jacketed Holder for S-1200 Series</p> <p>S-1200-111 Black Block To Set 0%T</p> <p>S-1200-512 Dust Cover for S-1200 Spectrophotometer</p> <p>S-1205-515 Halogen Bulb with Vertical Filament G4 12v20w for S-1205</p> <p>S-1200-505 Halogen Bulb, Package Of 1 (6v/10w, G4 Pin)</p>
<p>SPECIFICATION</p>	<p>Spectral Bandwidth 5 nm</p> <p>Monochromator Single Beam, Grating System 1200 lines/mm</p> <p>Wavelength Range 325-1000 nm</p> <p>Wavelength Accuracy</p>
<p>SIZE AND WEIGHT</p>	<p>Size 16 x 12 x 8 inch</p>



	(406 x 305 x 203mm)
Shipping Size	19.5 x 16.1 x 11.3 inch
	(496 x 410 x 286 mm)
Weight	18 lbs (8.2 kg)
Shipping Weight	21 lbs (9.5 kg)

