



Biohazard Shield

Model: 11410



Help prevent exposure to you and your staff from potentially harmful materials by splash or spill. Constructed of heavy 3/16" acrylic; our shields provide low cost - high value addition to your lab or practice. Use when opening urines, blood samples, or while pipeting and transferring hazardous liquids.



ADDITIONAL INFORMATION

SIZE AND WEIGHT	Size	11.5 x7.5 x 15 inch
	Shipping Size	13 x 8 x 16 inch
	Weight	2lbs
	Shipping Weight	3lbs

