


Phlebotomy Tray

Model: 48710



- 1 SINGLE PLACE TUBE-CUBE FOR ALCOHOL DISPENSER (PICTURED, NOT INCLUDED) OR ICE CUP
- 1 FOUR PLACE TUBE-CUBE FOR 50ML TUBES
- 2 NINE PLACE TUBE-CUBES FOR REGULAR BLOOD TUBES
- 5 BIN CUPS AND 2 RED SQUISHABLE HEARTS
- SHARPS CONTAINER AND TUBES NOT INCLUDED



The Flexi-Phleb has large carrying capacity without bulk. The tray provides a handy pull out drawer that is great for storing needles, bandages and other items. The top level accepts removable TubeCUBEs that secure a variety of tubes, ice cups for temperature sensitive samples, or sharps containers. Back at the lab, merely exchange your collected samples in Tube-CUBE



ADDITIONAL INFORMATION

MATERIAL	Wood
TOP SELECTION	No
SIZE AND WEIGHT	Size 9

