


# POLYSTYRENE, SQUARE, 10MM PATHLENGTH, 3.5ML CAPACITY, VISIBLE ONLY, BOX OF 500 PIECES

Model: S-90-302P-500



- POLYSTYRENE CUVETTE
- SQUARE 10MM PATHLENGTH
- 3.5ML CAPACITY
- CAN BE USED IN VISIBLE RANGE ONLY
- BOX OF 500 PIECES



Box of 500 polystyrene square cuvettes, 10mm pathlength, and 3.5ml capacity. Economical solution for high-volume visible spectrophotometric applications, offering reliable disposable cuvettes.



**ADDITIONAL INFORMATION**

<b>MATERIAL</b>	Wood
<b>TOP SELECTION</b>	No

