

# GLASS, RECTANGULAR, 50 MM PATHLENGTH, 17.5 ML CAPACITY VISIBLE ONLY, EACH

Model: S-90-321G



- GLASS CUVETTE
- RECTANGULAR 50MM PATHLENGTH
- 17.5ML CAPACITY CAN BE USED IN VISIBLE RANGE ONLY
- SINGLE



Glass rectangular cuvette with a 50mm pathlength and 17.5ml capacity, ideal for visible spectrophotometry. Perfect for measurements requiring a longer pathlength for higher absorbance readings.

## ADDITIONAL INFORMATION

<b>MATERIAL</b>	Wood
<b>TOP SELECTION</b>	No

