

GLASS, RECTANGULAR, 100 MM PATHLENGTH, 35.0 ML CAPACITY, VISIBLE ONLY, EACH

Model: S-90-322G



- GLASS CUVETTE
- RECTANGULAR 100MM PATHLENGTH
- 35.0ML CAPACITY
- CAN BE USED IN VISIBLE RANGE ONLY
- SINGLE



Glass rectangular cuvette with a 100mm pathlength and 35ml capacity, designed for visible range spectrophotometry. Suitable for larger sample volumes and concentrated solutions, offering precise absorbance readings.

ADDITIONAL INFORMATION

MATERIAL	Wood
TOP SELECTION	No

