

GLASS, SQUARE, 10X3 MM (DIA.), 0.07 ML, VISIBLE ONLY, EACH

Model: S-90-345FG



- GLASS CUVETTE
- FLOWCELL SQUARE WITH 10X3 MM (DIA.)AND 10MM PATHLENGTH
- 0.07 ML CAPACITY
- CAN BE USED IN VISIBLE RANGE ONLY
- SINGLE



Glass square cuvette with a 10x3mm diameter and a 0.07ml capacity, designed for low-volume samples in visible light spectrophotometry. Ideal for precise measurements of concentrated solutions.

ADDITIONAL INFORMATION

MATERIAL	Wood
TOP SELECTION	No

