


QUARTZ, SQUARE, 10MM PATHLENGTH, 3.5ML CAPACITY, UV-VIS , EACH

Model: S-90-309-1Q



- QUARTZ CUVETTE
- SQUARE 10MM PATHLENGTH
- 3.5ML CAPACITY
- CAN BE USED IN UV-VIS RANGE
- SINGLE



Single quartz cuvette, square shape with a 10mm pathlength and 3.5ml capacity. Designed for UV-Vis applications, ensuring accurate results across a wide spectral range, including ultraviolet measurements.

ADDITIONAL INFORMATION

MATERIAL	Wood
TOP SELECTION	No

