

QUARTZ, RECTANGULAR, 1 MM PATHLENGTH, 0.35 ML CAPACITY UV-VIS, EACH

Model: S-90-331Q



- QUARTZ CUVETTE
- RECTANGULAR 1 MM PATHLENGTH
- 0.35 ML CAPACITY
- CAN BE USED IN UV-VIS RANGE
- SINGLE



Quartz rectangular cuvette with a 1 mm pathlength and 0.35ml capacity, designed for UV-Vis spectrophotometry. Perfect for small volume and low concentration applications across the UV and visible spectrum.

ADDITIONAL INFORMATION

MATERIAL	Wood
TOP SELECTION	No

