


# QUARTZ, 10X2, 0.5ML, UV-VIS , EACH

Model: S-90-354Q



- QUARTZ CUVETTE
- 10X2MM CHAMBER AND 10MM PATHLENGTH
- 0.5ML CAPACITY
- CAN BE USED IN UV-VIS RANGE
- SINGLE



Quartz rectangular cuvette with dimensions of 10x2mm and a 0.5ml capacity, designed for UV-Vis spectrophotometry. Provides excellent transparency for precise measurements in both UV and visible light.

## ADDITIONAL INFORMATION

<b>MATERIAL</b>	Wood
<b>TOP SELECTION</b>	No

