


QUARTZ, LID, 10X1, 0.25ML, UV-VIS, EACH

Model: S-90-357Q



- QUARTZ CUVETTE
- WITH LID 10X1MM CHAMBER AND 10MM PATHLENGTH
- 0.25ML CAPACITY
- CAN BE USED IN UV-VIS RANGE
- SINGLE



Quartz cuvette lid with dimensions of 10x1mm and a 0.25ml capacity, designed for UV-Vis applications. Provides a secure closure for small-volume quartz cuvettes during spectrophotometric analysis.

ADDITIONAL INFORMATION

MATERIAL	Wood
TOP SELECTION	No

