


# QUARTZ, STOPPER, 10X2X2.5MM, 0.050ML, UV-VIS, EACH

Model: S-90-359Q



- QUARTZ CUVETTE
- STOPPER
- 10X2X2.5MM CHAMBER AND 10MM PATHLENGTH
- 0.050ML CAPACITY
- CAN BE USED IN UV-VIS RANGE
- SINGLE



Quartz cuvette stopper with dimensions of 10x2x2.5mm and a 0.050ml capacity, designed for UV-Vis spectrophotometry. Ideal for small-volume cuvettes requiring precise closure.

## ADDITIONAL INFORMATION

<b>MATERIAL</b>	Wood
<b>TOP SELECTION</b>	No

