


High School Microscope

Model: M220FLDM

- DUAL MONO HEAD
- WF10X EYEPIECE
- 4X,10X,40X DIN ACHROMAT OBJECTIVES
- ATTACHED MECHANICAL STAGE WITH 0.65 CONDENSER
- FLUORESCENT ILLUMINATOR



M220 Series Designed for High School biology, these durable and reliable microscopes are built to last and are backed by our Limited Lifetime Warranty. Select from a wide range of models and features to meet your needs. Pick from traditional models with tungsten bulbs or cool-lite fluorescent bulbs.

A selection of heads are available; chose from models with an inclined monocular head to dual, or "Teaching" heads. All models are parfocal and parcentered. All stages are drilled and tapped to accept optional mechanical stages; while certain models are available with mechanical stages already installed.

Each microscope includes arm-mounted cord wrap, dust cover and instruction manual.



ADDITIONAL INFORMATION

ACCESSORIES	Miscellaneous Accessories M220-7007 AC ADAPTER FOR FREEDOM JR. SERIES M220-7003 BULB, FLUORESCENT, 5 WATT M220-5101 GRADUATED MECHANICAL STAGE, ATTACHABLE </td>	
SPECIFICATION	Head Dual Head Eyepiece 10x WF (2) Objectives 4x, 10x, 40x Condenser NA 0.65 Diaphragm Disc, 5-hole Stage Mechanical Illumination Fluorescent Bulb 7 Watt Focusing Coarse/Fine	
SIZE AND WEIGHT	Actual Size 8.5 in x 5 in x 11.5 in (216 mm x 127 mm x 292 mm) Shipping Size 12 in x 8 in x 16 in (305 mm x 203 mm x 406 mm) Weight 8 lbs (3.6 kg) Shipping Weight 9 lbs (4.1 kg)	

