

# M260 SERIES MICROSCOPE

Model: M260

BINOCULAR, WF10X, 4X, 10X, 40XR, 100XR  
DIN, MECH., ABBE NA 1.25 & IRIS, LED



The M260 Binocular Microscope offers WF10X eyepieces, 4X, 10X, 40XR, and 100XR DIN objectives, a mechanical stage, and an Abbe condenser with NA 1.25 and iris diaphragm. With LED illumination, it ensures high-quality, bright light for precise specimen observation, making it an excellent choice for professional and educational applications.



ADDITIONAL INFORMATION

TOP SELECTION	Yes	
SPECIFICATION	Head	Binocular
	Eyepiece	WF10X
	Objectives	Achromat 4X/10X/40X/100X(oil)
	Condenser	NA 1.25
	Diaphragm	Iris
	Stage	Mechanical
	Illumination	LED
	Bulb	3 Watt LED
	Focusing	Coaxial
SIZE AND WEIGHT	Size	12.5 in x 8.5 in x 16.5 in
	Weight	13 lbs

