

# M260 SERIES MICROSCOPE

Model: M260

BINOCULAR, WF10X, 4X, 10X, 40XR, 100XR  
DIN, MECH., ABBE NA 1.25 & IRIS, LED



The M260 Binocular Microscope offers WF10X eyepieces, 4X, 10X, 40XR, and 100XR DIN objectives, a mechanical stage, and an Abbe condenser with NA 1.25 and iris diaphragm. With LED illumination, it ensures high-quality, bright light for precise specimen observation, making it an excellent choice for professional and educational applications.



**ADDITIONAL INFORMATION**

<b>TOP SELECTION</b>	Yes
<b>SPECIFICATION</b>	<b>Head</b> Binocular <b>Eyepiece</b> WF10X <b>Objectives</b> Achromat 4X/10X/40X/100X(oil) <b>Condenser</b> NA 1.25 <b>Diaphragm</b> Iris <b>Stage</b> Mechanical <b>Illumination</b> LED <b>Bulb</b> 3 Watt LED <b>Focusing</b> Coaxial
<b>SIZE AND WEIGHT</b>	<b>Size</b> 12.5 in x 8.5 in x 16.5 in <b>Weight</b> 13 lbs

