

DIGITAL POWERSPIN™ FX PRO

Model: C80508



- POWERSPIN™ FX PRO CENTRIFUGE
- TÜV RHEINLAND CERTIFIED
- FIXED SPEED 1000-3500 RPM
- 8 PLACE FIXED ANGEL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10 ML TUBES OR 16X15 ML TUBES



The C80508 PowerSpin™ FX PRO Centrifuge delivers consistent, fixed-speed operation at 3,400 RPM. Its 8-place rotor handles up to 8 x 10ml tubes, and the 99-minute digital timer ensures precise control. Certified by TÜV Rheinland, this centrifuge meets rigorous safety and EMC standards for reliable laboratory use.



ADDITIONAL INFORMATION

| | |
|------------------------|--|
| TOP SELECTION | Yes |
| ACCESSORIES | <p>Accessories</p> <p>Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-75</p> <p>C800-02 Regular tube shield, pack of 2</p> <p>C800-03 Adapter for short tubes, pack of 10</p> <p>C800-04 Pediatric or microtube adapter, pack of 2</p> <p>C800-05 Tube Lift Adapters, Pk/12</p> <p>C800-12 Balance Tube Holdster</p> |
| SPECIFICATION | <p>Number of Tubes 8</p> <p>Tube Capacity 8-10 ml</p> <p>Blood Yes</p> <p>Specimens Fecal Not recommended</p> <p>Urine Not recommended</p> <p>Microhematocrit Not recommended</p> <p>Speed Fixed</p> <p>Rotor Fixed Angle</p> <p>Maximum RPM 3500</p> <p>Maximum RCF(g) 1320</p> |
| SIZE AND WEIGHT | <p>Size 10.2 in x 11.8 in x 7.9 in</p> <p>Shipping Size 13.4 in x 14.8 in x 10.6 in</p> <p>Weight 13 lbs</p> <p>Shipping Weight 15 lbs</p> |
| WARRANTY | 1 years |

