

# DIGITAL POWERSPIN™ PORTA-SPIN™ CENTRIFUGE

Model: C82608



- POWERSPIN™ PX CENTRIFUGE
- TÜV RHEINLAND CERTIFIED
- VARIABLE SPEED 1000-3500 RPM
- 8 PLACE FIXED ANGEL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10 ML TUBES



The PowerSpin™ Porta-Spin™ C82608 portable centrifuge is optimized for on-the-go applications. Running on 12VDC, it's perfect for fieldwork or locations lacking standard power outlets. The unit features a fixed 3,500 RPM, an 8-position rotor designed for 10mL tubes, and a 99-minute digital timer for precise control. Meeting high safety and electromagnetic compatibility standards, and certified by TÜV Rheinland, this centrifuge provides dependable performance in diverse environments.



ADDITIONAL INFORMATION

TOP SELECTION	Yes	
ACCESSORIES	Accessories	
	C800-75	Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster
SPECIFICATION	Number of Tubes	8
	Tube Capacity	8-10 ml
	Blood	Yes
	Specimens Fecal	Not recommended
	Urine	Not recommended
	Microhematocrit	Not recommended
	Speed	Fixed
	Rotor	Fixed Angle
	Maximum RPM	3500
	Maximum RCF(g)	1397
SIZE AND WEIGHT	Size	10.2 in x 11.8 in x 7.9 in
	Shipping Size	13.4 in x 14.8 in x 10.6 in
	Weight	13 lbs
	Shipping Weight	15 lbs
WARRANTY	1 years	

