

Digital Powerspin™ LX PRO

Model: C85608

- POWERSPIN™ LX PRO CENTRIFUGE
- VARIABLE SPEED 1000-40000 RPM
- 8 PLACE FIXED ANGEL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10 ML TUBES




Short Description

- POWERSPIN™ LX PRO CENTRIFUGE
- VARIABLE SPEED 1000-40000 RPM
- 8 PLACE FIXED ANGEL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10 ML TUBES

The C85608 PowerSpin™ LX Centrifuge offers digital variable speed control from 300 to 4,000 RPM, providing flexibility for different sample types. With an 8-place rotor and a 99-minute digital timer, it can accommodate up to 8 x 15ml tubes, making it ideal for a wide range of laboratory applications. This centrifuge meets top safety and electromagnetic compatibility standards for reliable operation in laboratory environments. It complies with UL 61010-1:2012, CAN/CSA-C22.2 NO.61010-1-12, and IEC 61010-2-020:2016 for safety. It also adheres to FCC Part 15 Subpart B: 2017 and IECS-001: 2006 for electromagnetic compatibility, ensuring safe, interference-free use in labs globally.





Item	Icon	Name	Function
1		Blood program key	Press to load the parameter settings for blood separation.
2		Urine program key	Press to load the parameter settings for urine separation.
3		Fecal program key	Press to load the parameter settings for fecal separation.
4		SET key	Press to enter spin parameters
5		Plus key	Parameter increase
6		Minus key	Parameter decrease
7		ENTER key	Confirm the parameter.
8		Start/Stop key	Press to start or stop the spinning cycle

ADDITIONAL INFORMATION

TOP SELECTION	No
SHIPPING WEIGHT	22 lb
SHIPPING SIZE	Shipping Length 20 Shipping Height 16 Shipping Width 18
ACCESSORIES	<p>Accessories</p> <p>Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-75</p> <p>C8000-02 Regular Tube Shield, White, pack of 2</p> <p>C8000-03-8 5 ml Adapter for short tubes, pack of 8</p>
SPECIFICATION	<p>Number of Tubes 8</p> <p>Tube Capacity 8-15 ml</p> <p>Blood Yes</p> <p>Specimens Fecal Yes</p> <p>Urine Yes</p> <p>Microhematocrit With Adapters</p> <p>Speed Variable</p>



	Rotor	Fixed Angle
	Maximum RPM	4000
	Maximum RCF(g)	1320
SIZE AND WEIGHT	Shipping Size	20 in x 16 in x 18 in
	Weight	21 lbs
	Shipping Weight	22 lbs
WARRANTY		1 years

