

Digital Powerspin™ LX PRO

Model: C85608

- POWERSPIN™ LX PRO CENTRIFUGE
- VARIABLE SPEED 1000-40000 RPM
- 8 PLACE FIXED ANGEL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10 ML TUBES




Short Description

- POWERSPIN™ LX PRO CENTRIFUGE
- VARIABLE SPEED 1000-40000 RPM
- 8 PLACE FIXED ANGEL ROTOR
- 99 MIN DIGITAL TIMER
- CAPACITY 8X10 ML TUBES

The C85608 PowerSpin™ LX Centrifuge offers digital variable speed control from 300 to 4,000 RPM, providing flexibility for different sample types. With an 8-place rotor and a 99-minute digital timer, it can accommodate up to 8 x 15ml tubes, making it ideal for a wide range of laboratory applications. This centrifuge meets top safety and electromagnetic compatibility standards for reliable operation in laboratory environments. It complies with UL 61010-1:2012, CAN/CSA-C22.2 NO.61010-1-12, and IEC 61010-2-020:2016 for safety. It also adheres to FCC Part 15 Subpart B: 2017 and IECS-001: 2006 for electromagnetic compatibility, ensuring safe, interference-free use in labs globally.





ADDITIONAL INFORMATION

TOP SELECTION	No	
SHIPPING WEIGHT	22 lb	
SHIPPING SIZE	Shipping Length 20 Shipping Height 16 Shipping Width 18	
ACCESSORIES	Accessories	
	C800-75	Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case
	C8000-02	Regular Tube Shield, White, pack of 2
	C8000-03-8	5 ml Adapter for short tubes, pack of 8
SPECIFICATION	Number of Tubes	8
	Tube Capacity	8-15 ml
	Blood	Yes
	Specimens Fecal	Yes
	Urine	Yes
	Microhematocrit	With Adapters
	Speed	Variable



	Rotor	Fixed Angle
	Maximum RPM	4000
	Maximum RCF(g)	1320
SIZE AND WEIGHT	Shipping Size	20 in x 16 in x 18 in
	Weight	21 lbs
	Shipping Weight	22 lbs
WARRANTY		1 years

