


PRO SERIES 10 PLACE TUBE ROCKER

Model: L-TTR150

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING (± 2 RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER



Short Description

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING (±2RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER

The L-TTR150 Tube Rocker/Mixer can accommodate 10 tubes or vials and includes a pause button for convenient operation. Operating at 110V with 25 rocking cycles per minute, this model is designed for gentle and efficient mixing of small sample volumes in laboratory settings. The entire Pro Series line (L-TTR150, L-TTR250, and L-TTR450) meets the following safety and electromagnetic compatibility standards: UL 61010-1: 2012 R4.16 (US), CAN/CSA C22.2 No.61010-1-12+G11+G12 (Canada), FCC Part 15 Subpart B: 2017 (US), and ICES-001:2006 (Canada). These certifications ensure compliance with safety and EMC regulations for reliable operation in any laboratory environment.

ADDITIONAL INFORMATION

TOP SELECTION	Yes	
SPECIFICATION	Capacity	10 tubes
	Speed	Variable
	Voltage	110V
WARRANTY	1 years	

