


# PRO SERIES 16 PLACE TUBE ROCKER

Model: L-TTR250

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING ( $\pm 2$ RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER



## Short Description

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING ( $\pm 2$ RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER

The L-TTR250 PRO SERIES Tube Rocker/Mixer holds up to 16 tubes or vials and features a pause button for easy control. With 25 rocking cycles per minute and 110V operation, this model is ideal for medium-sized sample mixing in laboratories requiring reliable and consistent agitation. The entire Pro Series line (L-TTR150, L-TTR250, and L-TTR450) meets the following safety and electromagnetic compatibility standards: UL 61010-1: 2012 R4.16 (US), CAN/CSA C22.2 No.61010-1-12+G11+G12 (Canada), FCC Part 15 Subpart B: 2017 (US), and ICES-001:2006 (Canada). These certifications ensure compliance with safety and EMC regulations for reliable operation in any laboratory environment.

## ADDITIONAL INFORMATION

<b>TOP SELECTION</b>	No
<b>SPECIFICATION</b>	<b>Capacity</b> 16 tubes <b>Speed</b> Variable <b>Voltage</b> 110V
<b>WARRANTY</b>	1 years

