


PRO SERIES 16 PLACE TUBE ROCKER

Model: L-TTR250

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING (± 2 RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER



Short Description

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING (±2RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER

The L-TTR250 PRO SERIES Tube Rocker/Mixer holds up to 16 tubes or vials and features a pause button for easy control. With 25 rocking cycles per minute and 110V operation, this model is ideal for medium-sized sample mixing in laboratories requiring reliable and consistent agitation. The entire Pro Series line (L-TTR150, L-TTR250, and L-TTR450) meets the following safety and electromagnetic compatibility standards: UL 61010-1: 2012 R4.16 (US), CAN/CSA C22.2 No.61010-1-12+G11+G12 (Canada), FCC Part 15 Subpart B: 2017 (US), and ICES-001:2006 (Canada). These certifications ensure compliance with safety and EMC regulations for reliable operation in any laboratory environment.

ADDITIONAL INFORMATION

TOP SELECTION	No	
SPECIFICATION	Capacity	16 tubes
	Speed	Variable
	Voltage	110V
WARRANTY	1 years	

