


# PRO SERIES 30 PLACE TUBE ROCKER

Model: L-TTR450

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING ( $\pm 2$ RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER



## Short Description

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING ( $\pm 2$ RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER

The L-TTR450 PRO SERIES Tube Rocker/Mixer handles up to 30 tubes at 25 rocking cycles per minute. With a pause button and 110V operation, this model is ideal for high-throughput mixing of larger sample volumes in laboratory applications. The entire Pro Series line (L-TTR150, L-TTR250, and L-TTR450) meets the following safety and electromagnetic compatibility standards: UL 61010-1: 2012 R4.16 (US), CAN/CSA C22.2 No.61010-1-12+G11+G12 (Canada), FCC Part 15 Subpart B: 2017 (US), and ICES-001:2006 (Canada). These certifications ensure compliance with safety and EMC regulations for reliable operation in any laboratory environment.

## ADDITIONAL INFORMATION

|                      |  |
|----------------------|--|
| <b>TOP SELECTION</b> | No   |
| <b>SPECIFICATION</b> | <b>Capacity</b> 30 tubes<br><b>Speed</b> Variable<br><b>Voltage</b> 110V |
| <b>WARRANTY</b>      | 1 years  |

