


# PRO SERIES 30 PLACE TUBE ROCKER

Model: L-TTR450

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING ( $\pm 2$ RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER



Short Description

- SELECT FROM 32, 16 OR 10 TUBES MODELS
- DESIGNED FOR VARIOUS TUBE SIZES UP TO 16X120MM
- VARIABLE SPEED 10-30RPM, FOR GENTLE MIXING (±2RPM)
- NON-SLIP, REVERSIBLE AND AUTOCLAVABLE TUBE PAD, SILICON RUBBER

The L-TTR450 PRO SERIES Tube Rocker/Mixer handles up to 30 tubes at 25 rocking cycles per minute. With a pause button and 110V operation, this model is ideal for high-throughput mixing of larger sample volumes in laboratory applications. The entire Pro Series line (L-TTR150, L-TTR250, and L-TTR450) meets the following safety and electromagnetic compatibility standards: UL 61010-1: 2012 R4.16 (US), CAN/CSA C22.2 No.61010-1-12+GI1+GI2 (Canada), FCC Part 15 Subpart B: 2017 (US), and ICES-001:2006 (Canada). These certifications ensure compliance with safety and EMC regulations for reliable operation in any laboratory environment.

ADDITIONAL INFORMATION

TOP SELECTION	No	
SPECIFICATION	Capacity	30 tubes
	Speed	Variable
	Voltage	110V
WARRANTY	1 years	

