


Digital Powerspin™ LX

Model: C856D

- DIGITAL POWERSPIN™ LX CENTRIFUGE
- VARIABLE SPEED 300-4000 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. DIGITAL TIMER
- 6X10ML OR 3X15ML CAPACITY
- POWER INPUT 110V



Short Description

- DIGITAL POWERSPIN™ LX CENTRIFUGE
- VARIABLE SPEED 300-4000 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. DIGITAL TIMER
- 6X10ML OR 3X15ML CAPACITY
- POWER INPUT 110V

UNICO® Digital PowerSpin™ LX Centrifuge is a robust and user-friendly instrument designed for efficient sample separation. Its 6-place angled rotor accepts standard test tubes, making it suitable for a variety of laboratory procedures. With variable linear speed control up to 4000 RPM and a clear digital display, users can easily set and monitor the centrifugation process. Its compact design saves valuable bench space. Operates on 115V power.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories</p> <p>C800-00 8 Place Rotor</p> <p>C800-01 8 Place Rotor</p> <p>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p>



	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
	C800-12	Balance Tube Holdster
SPECIFICATION	Number of Tubes	6
	Tube Capacity	2-15 ml
	Blood	Yes
	Specimens Fecal	Yes
	Urine	Yes
	Microhematocrit	With Adapters
	Speed	Variable
	Rotor	Fixed
	Maximum RPM	300-4,000
	Maximum RCF(g)	Up to 1,828
SIZE AND WEIGHT	Size	10.2 in x 11.8 in x 7.9 in
	Shipping Size	13.4 in x 14.8 in x 10.6 in
	Weight	16 lbs
	Shipping Weight	18 lbs

