

CMH30 MOTOR

Model: CMH30-27



- INTEGRATED SPEED SENSOR: ENSURES ACCURATE RPM FEEDBACK AND PRECISE CONTROL.
- LARGE, DURABLE BRUSHES: EXTENDS MOTOR LIFESPAN AND MINIMIZES MAINTENANCE.
- MAINTAINS CENTRIFUGE PERFORMANCE: RESTORES OPTIMAL FUNCTIONALITY AND RELIABILITY.



The CMH30-27 Motor is a direct replacement for CMH30 microhematocrit centrifuges equipped with a digital timer. This 110V motor includes an integrated speed sensor for accurate RPM feedback and large, durable brushes for prolonged operational life. Restore your centrifuge's performance with this high-quality replacement part.



ADDITIONAL INFORMATION

| | |
|---------------|---------|
| TOP SELECTION | No |
| WARRANTY | 1 years |

