

# Infinity Microscope

Model: IP750

- SEIDENTOPF BINOCULAR HEAD
- WF10X PAIRED EYEPIECES
- 4X,10X,40XR,100XR DIN INFINITY ACHROMAT OBJECTIVES
- BUILT-IN MECHANICAL STAGE WITH ABBE NA1.25 CONDENSER & IRIS
- HALOGEN ILLUMINATION



The UNICO IP750 Infinity Binocular Microscope is a state-of-the-art instrument designed to deliver exceptional optical performance and precision for a range of applications. Ideal for educational, research, and clinical environments, this microscope combines advanced features to provide clear and accurate imaging.

Featuring WF10X widefield eyepieces, the IP750 offers a comfortable and expansive field of view, enhancing detailed observation and reducing eye strain during prolonged use. The microscope is equipped with a set of infinity achromatic objectives (4X, 10X, 40X, and 100X oil immersion), which are engineered to deliver superior image clarity and flat-field viewing across various magnifications. These objectives ensure minimal chromatic aberration and accurate specimen analysis.

The IP750 includes a 6V20W halogen illuminator with adjustable intensity control, allowing users to fine-tune the lighting conditions to suit different specimen types and enhance contrast. This feature ensures consistent and customizable illumination for optimal viewing.

With its high-quality optical components, precise engineering, and user-friendly design, the UNICO IP750 Infinity Binocular Microscope is an excellent choice for professionals seeking reliable performance and exceptional imaging quality for a wide range of microscopic applications.



## ADDITIONAL INFORMATION

TOP SELECTION	No
<b>ACCESSORIES</b>	<b>Bodies</b>
	IP730-3004 SEIDENTOPF BINOCULAR BODY WITHOUT EYEPIECES FOR IP730 SERIES MICROSCOPE
	IP730-3005 SEIDENTOPF TRINOCULAR BODY, 30 DEGREE INCLINED WITHOUT EYEPIECES FOR IP730 SERIES MICROSCOPE
	IP750-3004 SEIDENTOPF BINOCULAR BODY WITHOUT EYEPIECES FOR IP50 SERIES MICROSCOPE
	IP750-3005 SEIDENTOPF TRINOCULAR BODY, 30 DEGREE INCLINED WITHOUT EYEPIECES FOR IP750 SERIES MICROSCOPE
	<b>Optics</b>
	IP750-1502 10X EXTRA WIDEFIELD EYEPIECE FN 22MM FOR IP730 AND IP750 MICROSCOPE
	IP730-2101 4X DIN ACHROMAT OBJECTIVE, N.A. 0.10
	IP730-2102 10X DIN ACHROMAT OBJECTIVE, N.A. 0.25
	IP730-2103 40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS
	IP730-2104 100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, RETRACTABLE FRONT LENS
	IP730-8030 PHASE SET W/ PLAN PHASE 10X-40X-100X OBJECTIVES, TURRET CONDENSER WITH "0" POSITION
	IP750-2101 4X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2102 10X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2103 40X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2104 100X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2301 4X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE



	IP750-2302	10X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2303	40X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2304	100X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2422	INFINITY PLAN PHASE OBJECTIVE 10X FOR IP750 SERIES
	IP750-2423	INFINITY PLAN PHASE OBJECTIVE 40X FOR IP750 SERIES
	IP750-2424	INFINITY PLAN PHASE OBJECTIVE 100X FOR IP750 SERIES
	IP750-8030	PHASE SET W/ INFINITY PLAN PHASE 10X-40X-100X OBJECTIVES, TURRET CONDENSER WITH "0" POSITION
	<b>Miscellaneous Accessories</b>	
	IP750-5001	MECHANICAL STAGE WITHOUT ABBE CONDENSER FOR IP730 AND IP750 MICROSCOPE
	IP750-6002	PHASE CONTRAST CONDENSER, TURRET, NA1.25 WITH 10X-40X-100X PHASE RING AND "0" BRIGHT FIELD POSITION
	IP750-6003	DARKFIELD CONDENSER FOR IP730 AND IP750 SERIES MICROSCOPE
	IP730-8002	DUST COVER FOR IP730 AND 750 SERIES MICROSCOPE
	IP750-7006	6V 20W HALOGEN BULB, 4MM PIN FOR IP730 AND IP750 MICROSCOPE
	IP750-8901	HARD CARRYING CASE FOR IP730 AND IP750 SERIES MICROSCOPE
<b>SPECIFICATION</b>	<b>Head</b>	Binocular
	<b>Eyepiece</b>	WF10X
	<b>Objectives</b>	4x 10x 40XR 100XR Infinity Achr
	<b>Stage</b>	MECH.
	<b>Condenser</b>	Abbe N.A.1.25
	<b>Illuminator</b>	6V/20W /Halogen
<b>SIZE AND WEIGHT</b>	<b>Size</b>	10 in x 8.5 in x 15.4 in
	<b>Shipping Size</b>	11 in x 14 in x 18 in
	<b>Weight</b>	18 lbs



**Shipping Weight**

20 lbs

