

Infinity Microscope

Model: IP750

- SEIDENTOPF BINOCULAR HEAD
- WF10X PAIRED EYEPIECES
- 4X, 10X, 40XR, 100XR DIN INFINITY ACHROMAT OBJECTIVES
- BUILT-IN MECHANICAL STAGE WITH ABBE NA1.25 CONDENSER & IRIS
- HALOGEN ILLUMINATION



The UNICO IP750 Infinity Binocular Microscope is a state-of-the-art instrument designed to deliver exceptional optical performance and precision for a range of applications. Ideal for educational, research, and clinical environments, this microscope combines advanced features to provide clear and accurate imaging.

Featuring WF10X widefield eyepieces, the IP750 offers a comfortable and expansive field of view, enhancing detailed observation and reducing eye strain during prolonged use. The microscope is equipped with a set of infinity achromatic objectives (4X, 10X, 40X, and 100X oil immersion), which are engineered to deliver superior image clarity and flat-field viewing across various magnifications. These objectives ensure minimal chromatic aberration and accurate specimen analysis.

The IP750 includes a 6V20W halogen illuminator with adjustable intensity control, allowing users to fine-tune the lighting conditions to suit different specimen types and enhance contrast. This feature ensures consistent and customizable illumination for optimal viewing.

With its high-quality optical components, precise engineering, and user-friendly design, the UNICO IP750 Infinity Binocular Microscope is an excellent choice for professionals seeking reliable performance and exceptional imaging quality for a wide range of microscopic applications.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	Bodies
	SEIDENTOPF BINOCULAR BODY WITHOUT EYEPIECES FOR IP730 SERIES MICROSCOPE
	IP730-3004
	SEIDENTOPF TRINOCULAR BODY, 30 DEGREE INCLINED WITHOUT EYEPIECES FOR IP730 SERIES MICROSCOPE
	IP730-3005
	SEIDENTOPF BINOCULAR BODY WITHOUT EYEPIECES FOR IP50 SERIES MICROSCOPE
	IP750-3004
	SEIDENTOPF TRINOCULAR BODY, 30 DEGREE INCLINED WITHOUT EYEPIECES FOR IP750 SERIES MICROSCOPE
	IP750-3005
	Optics
	10X EXTRA WIDEFIELD EYEPIECE FN 22MM FOR IP730 AND IP750 MICROSCOPE
	IP750-1502
	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10
	IP730-2101
	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25
	IP730-2102
	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS
	IP730-2103
	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, RETRACTABLE FRONT LENS
	IP730-2104
	PHASE SET W/ PLAN PHASE 10X-40X-100X OBJECTIVES, TURRET CONDENSER WITH "0" POSITION
	IP730-8030
	4X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2101
	10X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2102
	40X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2103
	100X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2104
	4X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2301



	IP750-2302	10X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2303	40X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2304	100X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
	IP750-2422	INFINITY PLAN PHASE OBJECTIVE 10X FOR IP750 SERIES
	IP750-2423	INFINITY PLAN PHASE OBJECTIVE 40X FOR IP750 SERIES
	IP750-2424	INFINITY PLAN PHASE OBJECTIVE 100X FOR IP750 SERIES
	IP750-8030	PHASE SET W/ INFINITY PLAN PHASE 10X-40X-100X OBJECTIVES, TURRET CONDENSER WITH "0" POSITION
	Miscellaneous Accessories	
	IP750-5001	MECHANICAL STAGE WITHOUT ABBE CONDENSER FOR IP730 AND IP750 MICROSCOPE
	IP750-6002	PHASE CONTRAST CONDENSER, TURRET, NA1.25 WITH 10X-40X-100X PHASE RING AND "0" BRIGHT FIELD POSITION
SPECIFICATION	IP750-6003	DARKFIELD CONDERNSER FOR IP730 AND IP750 SERIES MICROSCOPE
	IP730-8002	DUST COVER FOR IP730 AND 750 SERIES MICROSCOPE
	IP750-7006	6V 20W HALOGEN BULB, 4MM PIN FOR IP730 AND IP750 MICROSCOPE
	IP750-8901	HARD CARRYING CASE FOR IP730 AND IP750 SERIES MICROSCOPE
SIZE AND WEIGHT	Head	Binocular
	Eyepiece	WF10X
	Objectives	4x 10x 40XR 100XR Infinity Achr
	Stage	MECH.
	Condenser	Abbe N.A.1.25
	Illuminator	6V/20W /Halogen
	Size	10 in x 8.5 in x 15.4 in
	Shipping Size	11 in x 14 in x 18 in
	Weight	18 lbs



	Shipping Weight	20 lbs
--	------------------------	--------

