


College Microscope

Model: M250



- SEIDENTOPF BINOCULAR HEAD
- WF10X PAIRED EYEPIECES
- 4X, 10X, 40XR, 100XR DIN ACHROMAT OBJECTIVES
- BUILT-IN MECHANICAL STAGE WITH ABBE NA1.25 CONDENSER & IRIS
- LED ILLUMINATION



The UNICO M250 Series microscopes are designed for college-level instruction and advanced discovery. This versatile series now features exclusively LED illumination, providing brighter, more energy-efficient lighting for optimal clarity and performance. Recent updates include a reverse nosepiece, a larger base, and a sturdier body construction, enhancing durability and value.

Choose from monocular, dual-view, or Seidentopf binocular head options, available in both corded and cordless models. All heads rotate 360° and feature diopter control on one eyepiece and adjustable interpupillary distance (55-75mm) for comfortable, personalized viewing.



ADDITIONAL INFORMATION

EQUIPMENT TYPE	Rotor, Adapter, Fuse
TOP SELECTION	No
ACCESSORIES	<p>Bodies</p> <p>M250-3001 MONOCULAR HEAD WITHOUT EYEPIECE</p> <p>M250-3002 DUAL VIEWING HEAD WITH DIOPTER WITHOUT EYEPIECES</p> <p>M250-3004 SEIDENTOPF BINOCULAR HEAD FOR M250 AND M250LED SERIES</p> <p>Optics</p> <p>M250-1401 10X WIDEFIELD EYEPIECE, SINGLE</p> <p>M250-1411 10X WIDEFIELD EYEPIECE WITH POINTER, SINGLE</p> <p>M250-1499 5X EYEPIECE, SINGLE</p> <p>M250-2101 4X DIN ACHROMAT OBJECTIVE, N.A. 0.10</p> <p>M250-2102 10X DIN ACHROMAT OBJECTIVE, N.A. 0.25</p> <p>M250-2103 40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS</p> <p>M250-2104 100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS</p> <p>Miscellaneous Accessories</p> <p>M250-5003 MECHANICAL STAGE WITH ABBE N.A. 1.25 CONDENSER AND IRIS DIAPHRAGM</p> <p>M250-7001 6V20W HALOGEN BULB, 4MM PIN FOR M250 SERIES</p> <p>M250-7003 FUSE 250V 1.5A</p> <p>M250-8001 DUST COVER FOR M250 SERIES</p> <p>M250-8902 EYEPIECE POINTER FOR M250 SERIES</p>
DOWNLOAD	M250.pdf



SPECIFICATION	Head	Binocular
	Eyepiece	10x WF (2)
	Objectives	4, 10, 40, 100x
	Condenser	NA1.25
	Diaphragm	Iris
	Stage	Mechanical
	Illumination	Halogen
	Focusing	Coaxial
SIZE AND WEIGHT	Actual Size	8.5 in x 6.5 in x 14.4 in
	Shipping Size	12.5 in x 8.5 in x 16.5 in
	Weight	13 lbs (6 kg)
	Shipping Weight	13 lbs (6 kg)
WARRANTY		1 years

