


Zoom Stereo Microscope

Model: ZM194

- BINOCULAR ZOOM STEREO HEAD
- WF10X PAIRED EYEPIECES
- 0.7X-4.5X CONTINUOUS ZOOM BODY
- FOCUSING STAND WITH TOP AND BOTTOM LED ILLUMINATION



The UNICO ZM194 Binocular Zoom Stereo Microscope is designed to meet the demands of detailed inspection, research, and educational applications with precision and ease. Combining advanced optics with a user-friendly design, this microscope provides exceptional clarity and versatility for a variety of microscopy needs.

Featuring 10X widefield eyepieces, the ZM194 offers a broad and clear field of view, reducing eye strain and enhancing user comfort during prolonged use. The microscope's zoom body provides a magnification range of 0.7X to 4.5X, allowing for smooth and continuous adjustment to achieve the optimal level of detail for any specimen.

The ZM194 is equipped with a sturdy rack stand and a large base, ensuring stability and support during extended observation sessions. The built-in top and bottom LED illumination systems provide bright, adjustable lighting, enhancing specimen visibility and contrast. These illumination features allow for customized lighting conditions to suit different specimen types and enhance overall image quality.

Whether used for industrial inspections, biological research, or educational purposes, the UNICO ZM194 Binocular Zoom Stereo Microscope offers reliable performance and high-quality imaging. Its combination of versatile magnification, stable design, and advanced LED lighting makes it an invaluable tool for a wide range of microscopy applications.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	Bodies
	D4-3001 Binocular Zoom Body 0.7X~4.5X
	D4-3002 Trinocular Zoom Body 0.7X~4.5X
	Eyepieces and Auxiliary Lens
	D4-1001 Widefield Eyepiece 10X FN: 19 mm
	D4-1002 Widefield Eyepiece 15X FN: 14.5 mm
	D4-1003 Widefield Eyepiece 20X FN: 11 mm
	D4-1004 Widefield Eyepiece 25X FN: 9 mm
	D4-1005 Extra Widefield Eyepiece 10X FN: 23 mm
	D4-2001 Auxiliary Lens 2X
	D4-2002 Auxiliary Lens 1.5X
	D4-2003 Auxiliary Lens 0.75X
	D4-2004 Auxiliary Lens 0.5X
	D4-2005 Auxiliary Lens 0.3X
	Stands, Arms and Bases
	D4-4001 PE Stand, pole type with extension pillar
	D4-4001HF PE-HF Stand, incident halogen (6V, 20W) and transmitted fluorescent illumination (5W)
	D4-4002 DPT Stand, with incident illuminator port
	D4-4003FF DFF Stand, incident and transmitted fluorescent illumination (5W)
	D4-4004 DE Stand, two-step extensible arm
	D4-4004HF Same as D4-4004 but with incident halogen (6V, 20W) and transmitted fluorescent illumination (5W)
D4-4006 Boom Stand with universal arm (D4-4101 included)	



	D4-4007	Boom Stand with 25 mm post carrier (D4-4106 included) and focusing mount (D4-4102 included)
	D4-4101	Universal arm (Bonder Arm) with 5/8 inch knuckle pin
	D4-4102	Focusing mount with 25 mm diameter through hole
	D4-4103	Probe arm
	D4-4104	Boom Stand with 5/8 inch knuckle hole
	D4-4105	Boom Stand with 25 mm post carrier (D4-4106 included)
	D4-4106	25 mm post carrier
	D4-4201	Diascopic base (only use with DPT stand D4-4002)
	Other Components and Accessories	
	D4-7002	Fiber optics dual-chrome gooseneck
	D4-7003	Fiber optics ring light
	D4-7004	Fiber optics light source 115V, 100W
	D4-8010	35 mm camera adapter and photo eyepiece system
	D4-8020	video adapter with external C thread, adjustable
	D4-8901	Objective protection ring with glass
	D4-8902	Reticule for 10X eyepiece (D4-1001)
	D4-8903	Reticule for WFH10X eyepiece (D4-1005)
SPECIFICATION	Head	Binocular
	Eyepiece	WF10X
	Body	0.7X~4.5X
	Stand	DF-HF Stand
	Illumination	Top halogen & Bottom fluorescent
SIZE AND WEIGHT	Actual Size	9.7" x 7" x 16.5"
	Shipping Size	10.8" x 14" x 17.5"
	Net Weight	15.5 lbs (7.0kg)
	Shipping Weight	17.5 lbs (7.7 kg)

