


Zoom Stereo Microscope

Model: ZM194



- BINOCULAR ZOOM STEREO HEAD
- WF10X PAIRED EYEPIECES
- 0.7X-4.5X CONTINUOUS ZOOM BODY
- FOCUSING STAND WITH TOP AND BOTTOM LED ILLUMINATION



The ZM191LED Binocular Zoom Stereo Microscope is a versatile and high-performance instrument ideal for a wide range of applications. It features 10x widefield eyepieces for comfortable viewing and a 0.7x-4.5x zoom body, providing versatile magnification options to suit your needs. The robust rack stand with a large base ensures stability during use, while the built-in top and bottom LED illumination offers bright and even lighting for clear observation of specimens.

**LIFETIME
LIMITED WARRANTY**



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	Bodies
	D4-3001 Binocular Zoom Body 0.7X~4.5X
	D4-3002 Trinocular Zoom Body 0.7X~4.5X
	Eyepieces and Auxiliary Lens
	D4-1001 Widefield Eyepiece 10X FN: 19 mm
	D4-1002 Widefield Eyepiece 15X FN: 14.5 mm
	D4-1003 Widefield Eyepiece 20X FN: 11 mm
	D4-1004 Widefield Eyepiece 25X FN: 9 mm
	D4-1005 Extra Widefield Eyepiece 10X FN: 23 mm
	D4-2001 Auxiliary Lens 2X
	D4-2002 Auxiliary Lens 1.5X
	D4-2003 Auxiliary Lens 0.75X
	D4-2004 Auxiliary Lens 0.5X
	D4-2005 Auxiliary Lens 0.3X
	Stands, Arms and Bases
	D4-4001 PE Stand, pole type with extension pillar
	D4-4001HF PE-HF Stand, incident halogen (6V, 20W) and transmitted fluorescent illumination (5W)
	D4-4002 DPT Stand, with incident illuminator port
	D4-4003FF DFF Stand, incident and transmitted fluorescent illumination (5W)
	D4-4004 DE Stand, two-step extensible arm
	D4-4004HF Same as D4-4004 but with incident halogen (6V, 20W) and transmitted fluorescent illumination (5W)
	D4-4006 Boom Stand with universal arm (D4-4101 included)



	D4-4007	Boom Stand with 25 mm post carrier (D4-4106 included) and focusing mount (D4-4102 included)
	D4-4101	Universal arm (Bonder Arm) with 5/8 inch knuckle pin
	D4-4102	Focusing mount with 25 mm diameter through hole
	D4-4103	Probe arm
	D4-4104	Boom Stand with 5/8 inch knuckle hole
	D4-4105	Boom Stand with 25 mm post carrier (D4-4106 included)
	D4-4106	25 mm post carrier
	D4-4201	Diascopic base (only use with DPT stand D4-4002)
	Other Components and Accessories	
	D4-7002	Fiber optics dual-chrome gooseneck
	D4-7003	Fiber optics ring light
	D4-7004	Fiber optics light source 115V, 100W
	D4-8010	35 mm camera adapter and photo eyepiece system
	D4-8020	video adapter with external C thread, adjustable
	D4-8901	Objective protection ring with glass
	D4-8902	Reticule for 10X eyepiece (D4-1001)
	D4-8903	Reticule for WFH10X eyepiece (D4-1005)
SPECIFICATION	Head	Binocular
	Eyepiece	WF10X
	Body	0.7X~4.5X
	Stand	DF-HF Stand
	Illumination	Top halogen & Bottom fluorescent
SIZE AND WEIGHT	Actual Size	9.7" x 7" x 16.5"
	Shipping Size	10.8" x 14" x 17.5"
	Net Weight	15.5 lbs (7.0kg)
	Shipping Weight	17.5 lbs (7.7 kg)
WARRANTY		1 years

