


DigiScout II Microscope Camera

Model: B6-8284



- INCLUDED C-MOUNT ADAPTER
- USB2.0 CABLE, MOUSE, AND DC-12V POWER ADAPTER
- USB & TF CARD OUTPUT



This innovative camera comes with a built-in Android operating system, providing a user-friendly interface that simplifies image capture and control. With a 1/3" sensor size and a Full HD resolution of 1024x600, the DigiScout II provides crisp, vibrant imaging at 30 FPS, perfect for capturing fine details of your specimens.



ADDITIONAL INFORMATION

TOP SELECTION	No
WARRANTY	1 years

