



DigiScout Microscope Tablet Camera.

Model: B6-8286


- INCLUDED C-MOUNT ADAPTER
- USB2.0 CABLE, MOUSE, AND DC-12V POWER ADAPTER
- USB & TF CARD OUTPUT



Short Description

- INCLUDED C-MOUNT ADAPTER
- USB2.0 CABLE, MOUSE, AND DC-12V POWER ADAPTER
- USB & TF CARD OUTPUT

Enhance your microscope experience with the DigiScout Microscope Tablet Camera. This camera features a 9-inch Full HD monitor integrated with a high-definition camera, designed for seamless observation and documentation. With a 1/3" sensor size and a Full HD resolution of 1024x600, the DigiScout provides crisp, vibrant imaging at 30 FPS, perfect for capturing fine details of your specimens



ADDITIONAL INFORMATION

TOP SELECTION	No
---------------	----

