

PowerSpin™ FX Centrifuge

Model: C806



- POWERSPIN™ FX CENTRIFUGE
- FIXED SPEED 3400 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML OR 3X15ML CAPACITY
- POWER INPUT 110V



The POWERSPIN™ FX Centrifuge is a practical and efficient solution for your centrifugation needs. This 110V centrifuge features a 6-place rotor, accommodating a variety of sample tubes. Its fixed speed of 3400 RPM and 30-minute timer simplify operation, while its durable construction ensures reliable performance in your laboratory. By using adapters, you can also spin tubes ranging from 2 – 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins™ are pre-drilled/tapped to accept optional Tube Holdsters™



ADDITIONAL INFORMATION

TOP SELECTION	No	
ACCESSORIES	Accessories (Sold Separately) C800-008 Place Rotor C800-018 Place Rotor Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case C800-75 C800-02Regular tube shield, pack of 2 C800-03Adapter for short tubes, pack of 10 C800-04Pediatric or microtube adapter, pack of 2 C800-05Tube Lift Adapters, Pk/12 C800-12Balance Tube Holdster	
SPECIFICATION	Number of Tubes 6 Tube Capacity 2-15 ml Blood Yes Specimens Fecal Not recommended Urine Not recommended Microhematocrit Not recommended Speed Fixed Rotor Fixed Angle Maximum RPM 3400 Maximum RCF(g) 1320	
SIZE AND WEIGHT	Size 10.2 in x 11.8 in x 7.9 in Shipping Size 13.4 in x 14.8 in x 10.6 in Weight 13 lbs Shipping Weight 15 lbs	



WARRANTY	1 years
----------	---------

