

PowerSpin™ FX Centrifuge

Model: C806H

- POWERSPIN™ FX CENTRIFUGE
- FIXED SPEED 3
- 400 RPM
- 6 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML OR 3X15ML CAPACITY
- WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 110V



The UNICO C806H PowerSpin™ FX Centrifuge is a top-tier choice for demanding laboratory environments, offering reliable performance and versatility. This advanced centrifuge features a 6-place rotor, designed to accommodate standard 8-15 ml tubes, and can also handle 2 - 7 ml tubes with optional adapters. With a fixed speed of 3,400 RPM, it ensures efficient separation for a range of materials, including blood and other biological samples.

The C806H model comes with an 18-place Tube Holdster™ rack, providing easy access and organization for up to 18 tubes or adapters. This optional side-mounted rack enhances your lab workflow by keeping essential tubes within reach and maximizing efficiency. The centrifuge is built with high-impact plastic housing for durability and features a mechanical 30-minute timer with a bell signal, ensuring consistent results. A safety switch in the lid halts the motor when the lid is raised, and a power indicator light shows when the unit is operating.

Compact yet powerful, the C806H measures 10.2 inches in width, 11.8 inches in depth, and 8 inches in height, weighing 15 lbs. Ideal for clinical labs, physician's offices, and chemical or industrial labs, the UNICO C806H PowerSpin™ FX Centrifuge delivers dependable performance and longevity for all your centrifugation needs.



ADDITIONAL INFORMATION

TOP SELECTION	No
ACCESSORIES	<p>Accessories</p> <p>C800-00 8 Place Rotor</p> <p>C800-01 8 Place Rotor</p> <p>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</p> <p>C800-02 Regular tube shield, pack of 2</p> <p>C800-03 Adapter for short tubes, pack of 10</p> <p>C800-04 Pediatric or microtube adapter, pack of 2</p> <p>C800-05 Tube Lift Adapters, Pk/12</p> <p>C800-12 Balance Tube Holdster</p>
SPECIFICATION	<p>Number of Tubes 6</p> <p>Tube Capacity 2-15 ml</p> <p>Blood Yes</p> <p>Specimens Fecal Not recommended</p> <p>Urine Not recommended</p> <p>Microhematocrit Not recommended</p> <p>Speed Fixed</p> <p>Rotor Fixed Angle</p> <p>Maximum RPM 3400</p> <p>Maximum RCF(g) 1320</p>
SIZE AND WEIGHT	<p>Size 10.2 in x 11.8 in x 7.9 in</p> <p>Shipping Size 13.4 in x 14.8 in x 10.6 in</p> <p>Weight 13 lbs</p> <p>Shipping Weight 15 lbs</p>

