


# PowerSpin™ FX Centrifuge

Model: C808

- POWERSPIN™ FX CENTRIFUGE
- FIXED SPEED 3400 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 8X10ML CAPACITY
- POWER INPUT 110V



UNICO® PowerSpin™ Series centrifuges are compact, economical and quiet. We offer both fixed (FX) and variable (LX) speed models with fixed angle rotors, for 6 or 8 tubes to spin 8- 15 ml tubes. By using adapters, you can also spin tubes ranging from 2 - 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins are pre-drilled/tapped to accept optional Tube Holdsters.



## ADDITIONAL INFORMATION

<b>ACCESSORIES</b>	<b>Accessories</b> <ul style="list-style-type: none"> <li>C800-00 8 Place Rotor</li> <li>C800-01 8 Place Rotor</li> <li>C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</li> <li>C800-02 Regular tube shield, pack of 2</li> <li>C800-03 Adapter for short tubes, pack of 10</li> <li>C800-04 Pediatric or microtube adapter, pack of 2</li> <li>C800-05 Tube Lift Adapters, Pk/12</li> <li>C800-12 Balance Tube Holdster</li> </ul>
<b>SPECIFICATION</b>	<ul style="list-style-type: none"> <li><b>Number of Tubes</b> 6</li> <li><b>Tube Capacity</b> 2-15 ml</li> <li><b>Blood</b> Yes</li> <li><b>Specimens Fecal</b> Not recomended</li> <li><b>Urine</b> Not recomended</li> <li><b>Microhematocrit</b> Not recomended</li> <li><b>Speed</b> Fixed</li> <li><b>Rotor</b> Fixed Angle</li> <li><b>Maximum RPM</b> 3400</li> <li><b>Maximum RCF(g)</b> 1320</li> </ul>
<b>SIZE AND WEIGHT</b>	<ul style="list-style-type: none"> <li><b>Size</b> 10.2 in x 11.8 in x 7.9 in</li> <li><b>Shipping Size</b> 13.4 in x 14.8 in x 10.6 in</li> <li><b>Weight</b> 13 lbs</li> <li><b>Shipping Weight</b> 15 lbs</li> </ul>

