


PowerSpin™ FX Centrifuge

Model: C808H

- POWERSPIN™ FX CENTRIFUGE
- FIXED SPEED 3400 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML CAPACITY
- WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 110V



UNICO® PowerSpin™ Series centrifuges are compact, economical and quiet. We offer both fixed (FX) and variable (LX) speed models with fixed angle rotors, for 6 or 8 tubes to spin 8– 15 ml tubes. By using adapters, you can also spin tubes ranging from 2 – 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins are pre-drilled/tapped to accept optional Tube Holdsters.



ADDITIONAL INFORMATION

| | |
|-----------------|--|
| ACCESSORIES | <div>Accessories</div> <div><div>C800-00</div><div>8 Place Rotor</div></div> <div><div>C800-01</div><div>8 Place Rotor</div></div> <div><div>C800-75</div><div>Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case</div></div> <div><div>C800-02</div><div>Regular tube shield, pack of 2</div></div> <div><div>C800-03</div><div>Adapter for short tubes, pack of 10</div></div> <div><div>C800-04</div><div>Pediatric or microtube adapter, pack of 2</div></div> <div><div>C800-05</div><div>Tube Lift Adapters, Pk/12</div></div> <div><div>C800-12</div><div>Balance Tube Holdster</div></div> |
| SPECIFICATION | <div><div>Number of Tubes</div><div>6</div></div> <div><div>Tube Capacity</div><div>2-15 ml</div></div> <div><div>Blood</div><div>Yes</div></div> <div><div>Specimens Fecal</div><div>Not recomend</div></div> <div><div>Urine</div><div>Not recomend</div></div> <div><div>Microhematocrit</div><div>Not recomend</div></div> <div><div>Speed</div><div>Fixed</div></div> <div><div>Rotor</div><div>Fixed Angle</div></div> <div><div>Maximum RPM</div><div>3400</div></div> <div><div>Maximum RCF(g)</div><div>1320</div></div> |
| SIZE AND WEIGHT | <div><div>Size</div><div>10.2 in x 11.8 in x 7.9 in</div></div> <div><div>Shipping Size</div><div>13.4 in x 14.8 in x 10.6 in</div></div> <div><div>Weight</div><div>13 lbs</div></div> <div><div>Shipping Weight</div><div>15 lbs</div></div> |

