



PowerSpin™ FX Centrifuge

Product Images



Short Description

- PowerSpin™ FX CENTRIFUGE
- FIXED SPEED 3400 RPM
- 8 PLACE FIXED ANGLE ROTOR
- 30 MIN. MECHANICAL TIMER
- 6X10ML CAPACITY
- WITH 18 PLACE TUBE HOLDSTER™ RACK
- POWER INPUT 110V

Description

UNICO® PowerSpin™ Series centrifuges are compact, economical and quiet. We offer both fixed (FX) and variable (LX) speed models with fixed angle rotors, for 6 or 8 tubes to spin 8– 15 ml tubes. By using adapters, you can also spin tubes ranging from 2 – 7 ml. An optional Tube Holdster™ allows you to have up to 18 tubes or adapters at your fingertips, side-mounted on the centrifuge. PowerSpins are pre-drilled/tapped to accept optional Tube Holdsters.

Additional Information

Accessories	Header Accessories C800-00 8 Place Rotor C800-01 8 Place Rotor C800-75 Tachometer, Laser Digital with Reflective Tape, Batteries, In A Storage Case C800-02 Regular tube shield, pack of 2 C800-03 Adapter for short tubes, pack of 10 C800-04 Pediatric or microtube adapter, pack of 2 C800-05 Tube Lift Adapters, Pk/12 C800-12 Balance Tube Holdster
Specification	Number of Tubes 6 Tube Capacity 2-15 ml Blood Yes Specimens Fecal Not recommended Urine Not recommended Microhematocrit Not recommended Speed Fixed Rotor Fixed Angle Maximum RPM 3400 Maximum RCF(g) 1320
Size and Weight	Size 10.2 in x 11.8 in x 7.9 in Shipping Size 13.4 in x 14.8 in x 10.6 in Weight 13 lbs Shipping Weight 15 lbs

